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Handling:

Handle all glass filter standards using only the metal frames. Do not touch the glass surfaces. Any contamination of the surfaces can result in erroneous results in calibration and subsequent spectrophotometer usage.

Be consistent in filter orientation each time the standard is used. Measurement differences may sometimes be seen depending on which way the Standard is inserted into the instrument, as the filter is not at the center of the sample compartment. GEX recommends facing the glass portion into the main body of the instrument.

Refer to the product manual for more information.

A GEX Corporation Dosimetry Procedure Template is located on the GEX website www.gexcorp.com which provides a detailed procedure for the use of the standards along with a template Form for recording test results (GEX Document# 100-262). We suggest using this procedure and form in place of that provided by the manufacturer.

Calibration:

GEX recommends that the Standards be calibrated every two years, but the specific calibration interval should be defined by the end-user. However, the Standards must be recertified if you suspect that they are damaged or no longer exhibit the certified values. This service (catalog # S4320) is available through GEX. Please contact GEX Corporation for instructions before returning the Standards.

Warranty:

Please fill out and send the warranty card to Thermo Electron Corporation. GEX Corporation and Thermo Electron Corporation do not warrant user modifications or the use of unapproved cleaning methods.